



LORISIAN SILVER TEST

What will the test find out?

This test will analyse your IgG antibody reactions to 156 food ingredients, the laboratory will test your IgG antibody reactions to 156 food ingredients to identify your personal 'food fingerprint'.

The LorisiaN Silver Test will test your IgG antibody reaction to the foods on reverse of this page.

How does the test work?

The test is simple to take and can be used at home or taken at the club. All that is required is a small blood sample, which involves a simple prick of the finger with a small, sterilised lancet. Once the blood sample has been taken it is then posted straight to the laboratory for analysis.

How accurate are the results?

The LorisiaN laboratory is regularly audited by independent authorities and uses expert staff who are fully qualified, trained and experienced in using the latest scientific equipment to ensure the utmost accuracy of testing.

The LorisiaN laboratory has a >98% reproducibility rate. This means that, on average, two separate tests carried out with the same blood sample (taken on the same day) will correlate with each other >98% of the time.

What is the next step?

If you suffer from any of the following symptoms:

- IBS
- Bloating
- Digestive Problems
- Migraine
- Headaches
- Joint Pain
- Itchy Skin
- Eczema
- Fatigue
- Low Mood

You may benefit from taking a LorisiaN Food Intolerance Test.

To book in for your test, please speak to a member of staff in reception or email us on:

info@breedonprioryhealthclub.co.uk